Needleholders Nadelhalter

Needleholders

There is more to a good needleholder than meets the eye.

This is especially true for S&T's line:

- **Secure grip** on even the smallest needle and finest suture material is possible due to parallel closure over the entire jaw length.
- **Easy tying** of the finest suture material is possible due to perfectly designed jaws.
- Closing pressure is precisely calibrated to facilitate holding of the instrument while **minimizing hand fatigue**.
- **Safe suture handling** as all sharp corners and rough spots are eliminated.

Selection of the proper Needleholder

In addition to the observations described in the beginning pages of the catalog, there are four points that you should take into consideration when selecting a needleholder:

- your preference for the needleholder tip, straight or curved
- your preference for locking or non-locking needleholder
- your preference for a hinged or non-hinged needleholder
- Round or flat handle

Needleholder tip

Curved tip needleholders are the most commonly used instruments, however on occasion some surgeons prefer straight needleholders.

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
S&T®

S&T® Needleholders *Nadelhalter*

Needleholders with the designation "C" in the article code are based on the design by Castroviejo and are basically fine, flat handle needleholders. The designation "B" stands for round handle needleholders which were introduced by Barraquer.

Both designs originate from ophthalmic instrument designs.

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Prof. Dr. Ignacio Barraquer y Barraquer, Barcelona, Spanien

Needleholders with Lock

The lock on a needleholder is mainly intended to be a convenient means for parking the needle, and passing it to the surgeon. The lock is not intended to be an aid for holding the needle while suturing.

S&T needleholders with a lock are carefully designed and calibrated. The lock may be disengaged and re-engaged without the inconvenience of letting go of the needle. Engagement and disengagement is marked by an audible click. Since the finger pressure required to release the lock may cause unwanted movement at a critical moment. We recommend that the passing of the needle through tissue should be done with the lock not engaged.

Note

Refrain from using microsurgical needle holders to hold a macrosurgical needle. Needles over 8 mm length are so thick and the force required to grip and pass them is so high, that even a single use may permanently damage the jaws of a micro needleholder.

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


15 cm long, round handle 8 mm dia

Länge 15 cm, Rundgriff, 8 mm Ø


BS-15-8	—	Art. No. 00551	Ø 0.40	●				○	●	●	●	●	●
B-15-8	⤿	Art. No. 00552	Ø 0.40	●				○	●	●	●	●	●
BL-15-8	⤿	🔒	Art. No. 00553	Ø 0.40	●	●		○	●	●	●	●	●
B-15-8.3	⤿	Art. No. 00763	Ø 0.30	●				○	●	●	●	●	●

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**Non-Hinged
Needleholders**

(angulated forceps)

Angulated forceps are extremely useful light-duty needle holders, very widely preferred for lab work. They are also useful for fine clinical work on vessels up to 2 mm in diameter. Above this size a hinged needle holder is recommended, since a firmer grip on the needle is required when working with thicker tissue.

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Needleholder forceps

Needleholder forceps have stronger tips than standard microsurgical forceps. Needleholder forceps are produced to the same rigid specifications of our spring-handled needleholder. **S&T** needleholder forceps are designed to facilitate rotational technique.

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
Needleholder forceps, curved,
15 cm long, round handle 8 mm dia,
Model registered

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
Needleholder forceps, 15 cm long,
round handle, curved 45°, 8 mm dia,
Model registered

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Needleholder forceps, curved,
18 cm long, round handle 8 mm dia,
Model registered

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Nadelhalter

Double action needleholders

Double action needleholders are especially advantageous in situations where access to an operation site is limited. Double action needleholders feature a reduced handle width compared to a conventional spring handled needleholder.

These instruments are an original **S&T** development.

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