GEM

Global Excellence in Microsurgery

GEM Microvascular Anastomotic COUPLER Device and System

EN	Instructions for use.....	6
FR	Instructions d'utilisation.....	16
DE	Gebrauchsanweisung.....	27
IT	Istruzioni per l'uso.....	39
ES	Instrucciones de uso.....	50
NL	Gebruiksaanwijzing.....	61
DA	Betjeningsvejledning.....	73
SV	Användarinstruktioner.....	83
TR	Kullanım Talimatları.....	93
EL	Οδηγίες χρήσης.....	103
JA	使用説明書.....	115
ZH	使用说明书.....	125
NO	Instruksjoner for bruk.....	134
PT	Instruções de Uso.....	144
PL	Instrukcja użycia.....	155
KO	사용설명서.....	167

Microvascular Anastomotic COUPLER Device and System
Appareil et système Microvasculaire Anastomotique COUPLER
Microvascular Anastomotic COUPLER-Gerät und System
Dispositivo e sistema COUPLER microvascolare anastomotico
Sistema y dispositivo microvascular anastomótico COUPLER
Microvasculair anastomotisch COUPLER apparaat en systeem
Mikrovaskulær anastomotisk COUPLER enhed og system
Mikrovaskulär COUPLER-anordning och -system för anastomos
Mikrovasküler Anastomotik COUPLER Cihaz ve Sistemi
Σύστημα και συσκευή COUPLER μικροαγγειακής αναστόμωσης
微小血管吻合用COUPLERデバイスおよびシステム
微血管吻合 COUPLER 简体中文和系统
Mikrovaskulær anastomotisk COUPLER-enhet og system
Sistema e Dispositivo Anastomótico Microvascular COUPLER
Urządzenie i system COUPLER do wykonywania zespoń mikronaczyniowych
미세혈관 문합용 COUPLER 기기 및 시스템

CE
2797


Figure 1

Figure 2

Figure 3


Figure 4


Figure 5a


Figure 5b


Figure 6

Figure 7


Figure 8


Figure 9


Figure 10


Figure 11


Figure 12


Figure 13


Figure 14


Figure 15


Figure 16


Figure 17


Figure 18


Figure 19


Figure 20a


Figure 20b


Figure 21a


Figure 21b


Figure 22


Figure 23

SYMBOL GLOSSARY PER US FD&C ACT:

Standard	Symbol	Symbol Title	Symbol Meaning	Symbol Number
ISO-15223-1*		Do not reuse	Do not reuse	5.4.2
ISO-15223-1		Use by Date	Use by Date	5.1.4
ISO-15223-1	STERILE R	Sterilized using irradiation	Sterilize using irradiation	5.2.4
ISO-15223-1		Do not use if package is damaged	Do not use if the product sterilization barrier or its packaging is compromised.	5.2.8
ISO-15223-1		Consult instructions for use	Consult instructions for use	5.4.3
ISO-15223-1		Caution	Attention, consult accompanying documents.	5.1.5
			This product and package is not made with natural rubber latex.	

Standard	Symbol	Symbol Title	Symbol Meaning	Symbol Number
	CONTENT		Content	
	Rx Only		CAUTION: Federal (USA) Law restricts this device to sale by or on the order of a physician.	
ISO-15223-1	REF	Catalogue number	Catalogue number	5.1.6
ISO-15223-1	LOT	Batch code	Lot number	5.1.5
ISO-15223-1		Manufacturer	Manufacturer	5.1.1
ISO-15223-1	EC REP	Authorized representative in the European Community	Authorized representative in the European Community	5.1.2
ISO-15223-1	MD	Medical device	Indicates the item is a medical device	5.7.4
ISO-15223-1		Do not re-sterilize	Do not re-sterilize	5.2.6
ISO-15223-1		Date of manufacture	Date of manufacture	5.1.3

*ISO 15223-1:2016 Medical Devices – Symbols to be used with medical device labels, labeling and information to be supplied – Part 1: General requirements

Additional symbols and graphics on the product labeling that are not required by the US FD&C Act:

Symbol	Symbol Description
	Made in the U.S.A.
	Manufacturer part number
	Internal tracking number
	Serial number

Graphic	Graphic Meaning
	Size of the GEM COUPLER Device (inner diameter of the COUPLER rings)

DESCRIPTION:

The Synovis MCA GEM Microvascular Anastomotic COUPLER Device and System have been specifically designed for use in the anastomosis of blood vessels.

The GEM Microvascular Anastomotic COUPLER Device rings are made of polyethylene and surgical grade stainless steel pins. A protective cover and jaw assembly protect the rings and allow for easy loading onto the Anastomotic Instrument. Both the protective cover and jaw assembly are disposable. The GEM Microvascular Anastomotic COUPLER Device is single-use and available in various sizes.

The GEM Microvascular Anastomotic COUPLER System consists of a reusable Anastomotic Instrument (surgical-grade stainless steel and titanium), reusable COUPLER Forceps

(surgical-grade stainless steel), a reusable Double-Ended Vessel Measuring Gauge (surgical-grade stainless steel), and a Sterilization Tray (anodized aluminum). The GEM Microvascular Anastomotic COUPLER System is available as a complete system or as single items.

INDICATIONS FOR USE:

The GEM Microvascular Anastomotic COUPLER Device is intended to be used in the anastomosis of veins and arteries normally encountered in microsurgical procedures only in the peripheral vascular system.

The GEM Microvascular Anastomotic COUPLER Device is intended for use with veins and arteries having an outside diameter no smaller than 0.8 mm and no larger than 4.3 mm and a wall thickness of 0.5 mm or less. The degree of vessel spasm and the elasticity of the vessel should be considered when choosing the COUPLER Device size to be used. Use a GEM Vessel Measuring Gauge to approximate the OUTER diameter of the vessel for selection of an appropriate COUPLER Device size.

NOTE: Approximate the vessel to an appropriate COUPLER size using a GEM Vessel Measuring Gauge. In an end-to-end anastomosis, the two vessel ends should be approximately the same size as the inside diameter of the COUPLER Device being used. In an end-to-side anastomosis, the opening made to the side vessel should be approximately the same size as the inside diameter of the COUPLER Device being used. If the end vessel size is outside of the specified range, do not use the COUPLER Device to mate the vessels.

EN

CONTRAINDICATIONS:

The COUPLER Device is not indicated for use in patients presenting conditions that would normally preclude microvascular repair with suture technique. Examples of such conditions include, but are not limited to:

- Pre-existing or suspected peripheral vascular disease,
- Ongoing irradiation of the area of reconstruction,
- Clinical infection of the area of reconstruction,
- Anticipated infection due to significant contamination of the area of reconstruction,
- Friability of the vascular tissue due to sclerotic conditions,
- Concurrent diabetes mellitus, or
- Concurrent corticosteroid therapy

WARNINGS:

- Failure to use the Vessel Measuring Gauge to approximate the vessel size could result in using a COUPLER Device of an inappropriate size. Using a ring too large for the vessel may result in stressing or tearing of the vessel wall and a compromised anastomosis. Using a ring too small for the vessel may unduly constrict the vessel and lead to thrombosis or ring separation.
- Failure to squeeze the COUPLER jaws with a hemostat or similar instrument prior to ejection of the joined rings may result in an inadequate friction fit and possible ring separation. **Inspect the anastomotic site** to ensure that the anastomosis has been satisfactorily completed.

- The COUPLER Device is supplied sterile and is single use only. **Do not resterilize or reuse** the COUPLER Device.
 - Resterilization may compromise the structural integrity of the product which may lead to incomplete anastomosis.
 - Device cannot be reused due to possible structural damage during first use which may lead to incomplete anastomosis.
- Do not use the COUPLER Device if the package has been opened or appears to be damaged or compromised as sterility may be compromised. Failure to observe this warning may result in surgical infection.
- Safe use of the COUPLER Device for the anastomosis of tubular structures other than veins and arteries has not been established.
- Safe use of the COUPLER Device for the anastomosis of growing vessels in children or adolescents has not been established. Not intended for fetal use.
- Security of an anastomosis utilizing COUPLER Devices that have been approximated, reopened, and then reapproximated has not been demonstrated. When reapproximation of the anastomosis is desired, the vessel should be removed from each ring and a new COUPLER Device utilized.
- The Anastomotic Instrument, Vessel Measuring Gauge, COUPLER Forceps, and Sterilization Tray **must be sterilized prior to use.**
- The Anastomotic Instrument, Vessel Measuring Gauge, COUPLER Forceps, and Sterilization Tray should be

thoroughly inspected before use. Instruments that are damaged and/or in need of repair should not be used.

- When performing an end-to-side anastomosis with this COUPLER device, the lumen of the “side” vessel is narrowed slightly. For this reason, the diameter of the “side” vessel should be larger than that of the “end” vessel when completing such a procedure. The opening made to the side vessel should be approximately the same size as the inside diameter of the COUPLER device being used.

CAUTIONS:

- Use of the COUPLER Device involves potential risks normally associated with any implanted device, e.g., infection, perforation, or laceration of vessels, erosion, implant rejection, or device dislodgement/migration.

INSTRUCTIONS FOR USE:

These Instructions for Use are designed for proper use of this device. They are not intended to serve as a reference to surgical technique, to supersede institutional protocols or professional clinical judgment regarding patient care.

It is the responsibility of the clinician to inform the patient that he/she is the recipient of permanent implants which contain metal components (surgical-grade stainless steel pins). The COUPLER Device rings have been evaluated with a 1.5 Tesla and a 3.0 Tesla magnetic field and no change in displacement was observed in each of three orthogonal planes.^{1,3} The stainless steel pins in the COUPLER Devices are nominally nonferromagnetic. However, the US Food and Drug Administration (FDA) has made

recommendations for any medical device implanted which have metallic components to include:

- Documentation in the official medical record of the identity of the implant (manufacturer, model number, lot and serial numbers, and identifying marks, if any).
- Documentation of the technique and results of any magnetic testing performed on the implant or that no such testing was done.
- Patient education regarding the particular implant and recommendation for identifying medical alert card, bracelet, or necklace characterizing the implanted device.²

3.0mm COUPLER Size or Smaller:

End-to-End Anastomosis:

Using conventional microsurgical technique, mobilize a minimum of 1 cm of each vessel end. Using vascular clamps, clamp off the vessel(s) and irrigate the vessel openings. The COUPLER requires a greater amount of free vessel within the clamps than a conventional suture repair.

1. After gentle dilation, estimate the **outer** diameter of each vessel using the Vessel Measuring Gauge. The circular guides on the gauge **should not** be placed inside the vessel lumen (See Figure 1). If there is a size discrepancy between the two vessels, use the measurement of the smaller vessel to choose the appropriate COUPLER Device. The degree of vessel spasm and the elasticity of the vessel should be considered when choosing the COUPLER Device size to be used.

2. Select the appropriate size COUPLER Device. Both vessel ends should be approximately the same size as the inside diameter of the COUPLER Device being selected.
3. Remove the lid from the outer tray and aseptically remove the inner tray; the inner tray may be placed in the sterile field. Inspect the inner tray. Do not use if the inner tray is damaged or if the seals are not intact. Remove the lid from the inner tray.
4. Turn the Anastomotic Instrument knob fully counterclockwise, and then insert the COUPLER onto the Anastomotic Instrument. **The matching indicator arrows on the COUPLER Device and the Anastomotic Instrument should be pointing toward each other when loading (See Figures 2 & 3). Ensure that an audible click is heard for proper loading.**
5. Remove COUPLER Device from the protective cover by pulling it firmly away from the Anastomotic Instrument (See Figure 4).
6. Visually inspect to see that both rings are seated at the bottom of the U portion of the jaw and the pins are not bent (See Figures 5a & 5b). If pins are bent, do not attempt to straighten. Instead use a new COUPLER Device.
7. Place the Anastomotic Instrument perpendicular to the vessel(s), with the COUPLER jaw assembly near the two vessel ends. Pull one vessel end through one of the COUPLER rings using microsurgical forceps (See Figure 6).
8. Take a bite of approximately one to two pin diameters of the vessel wall and intimal lining, evert 90 degrees and impale

onto one pin. Proceeding in a triangular fashion, impale the vessel firmly upon every other pin, completing three pins (See Figure 7). Complete vessel placement on the ring by impaling the vessel upon the remaining three intermediate pins (See Figure 8). Ensure that both the vessel wall and the intimal layer are fully impaled upon each pin to reduce the risk of thrombosis. Should the vessel wall tear during impalement, remove the vessel, trim the end, and repeat the procedure. For examples of improper impalement of the vessel see Figure 9.

9. Repeat Steps 7 and 8 to impale the other vessel end upon the second COUPLER ring.
10. When both vessel ends have been suitably impaled, visually inspect to ensure that both rings are seated at the bottom of the U portion of the jaw and the pins are not bent (See Figures 5a and 5b). Bring the rings together (See Figures 10 & 11) by turning the Anastomotic Instrument knob clockwise. Turn the knob only until the ejector rod has just started to move the now joined rings.
11. **Prior to ejecting the joined rings, gently squeeze the end of the apposed jaws with a small hemostat (See Figure 12) to ensure ring approximation and a tight friction fit.** Turn the Anastomotic Instrument knob further clockwise to eject the joined rings.
12. Check the anastomosis under the operating microscope before opening the vascular clamps. Remove the clamps and **inspect the anastomotic site to ensure that the anastomosis has been satisfactorily completed** (patent vessel without leakage).

13. To remove the jaw assembly turn the Anastomotic Instrument knob fully counterclockwise (See Figure 13). Press the release button, located near the arrow on the Anastomotic Instrument, and remove the jaw assembly (See Figure 14).
14. Rinse the Anastomotic Instrument with water after use.

3.5mm COUPLER Size or Larger:

End-to-End Anastomosis:

1. to 7. Follow the same directions as for 3.0mm COUPLER Size or Smaller End-to-End Anastomosis (Steps 1 through 7).
8. Take a bite of approximately one to two pin diameters of the vessel wall and intimal lining, evert 90 degrees and impale onto the pin situated nearest to the open part of the Jaw Assembly (open end of the U portion of the jaw). Impale the opposite side of the vessel opening to the pin directly across from the initial pin. Next, impale the vessel onto the pins located near the sides of the ring, keeping the vessel as evenly spaced as possible between the four pins (See Figure 15). Continue vessel placement on the ring by impaling the vessel onto the two remaining pins near the open end of the Jaw Assembly. Complete by impaling the vessel onto the last two pins near the bottom of the Jaw Assembly (bottom of the U portion of the jaw); this final step prevents the ring from sliding out of the Jaw Assembly prematurely (See Figure 16). Ensure that both the vessel wall and the intimal lining are fully impaled upon each pin to reduce the risk of thrombosis. Should the vessel wall tear during impalement, remove the vessel, trim the end, and repeat the procedure. For examples of improper impalement of the vessel see Figure 17.

9. to 14. Follow the same directions as for 3.0mm COUPLER Size or Smaller End-to-End Anastomosis (Steps 9 through 14).

All COUPLER Sizes:

End-To-Side Anastomosis:

Using conventional microsurgical technique, mobilize a minimum of 1 cm of the “end” vessel. Clamp-off the vessel and irrigate the vessel lumen. Mobilize a minimum of 2 cm of the “side” vessel and clamp off the vessel.

1. When performing an end-to-side anastomosis with this COUPLER device, the lumen of the “side” vessel is narrowed slightly. For this reason, the diameter of the “side” vessel should be larger than that of the “end” vessel when completing such a procedure.
2. Estimate the outer diameter of the “end” vessel using the Vessel Measuring Gauge. The circular guides on the gauge **should not** be placed inside the vessel lumen (See Figure 1).
3. Select the appropriate size COUPLER Device.
4. Remove the lid from the outer tray and aseptically remove the inner tray; the inner tray may be placed in the sterile field. Inspect the inner tray. Do not use if the inner tray is damaged or if the seals are not intact. Remove the lid from the inner tray.
5. Turn the Anastomotic Instrument knob fully counterclockwise, and then insert the COUPLER onto the Anastomotic Instrument. **The matching indicator arrows on the COUPLER Device and the Anastomotic Instrument should be pointing**

EN

toward each other when loading (See Figures 2 & 3). **Ensure that an audible click is heard for proper loading.**

6. Remove COUPLER Device from the protective cover by pulling it firmly away from the Anastomotic Instrument (See Figure 4).
7. Visually inspect to see that both rings are seated at the bottom of the U portion of the jaw and the pins are not bent (See Figures 5a and 5b). If pins are bent, do not attempt to straighten. Instead use a new COUPLER Device.
8. Place the Anastomotic Instrument handle perpendicular to the direction of the “end” vessel. Place the “end” vessel upon one ring as described in Steps 7 and 8 for End-to-End Anastomosis directions of the appropriate COUPLER size.
9. Create a transverse incision in the “side” vessel of a length no greater than the internal diameter of the COUPLER device selected. Slide the clamps together slightly to remove tension and to open the incision (See Figure 18). Irrigate the vessel lumen through the created opening.
10. Using microsurgical forceps, grasp the vessel wall and intimal lining near one end of the transverse incision and pull them through the remaining ring. Evert the vessel wall and intimal lining 180 degrees and impale the vessel first upon the pins situated nearest the incision end (See Figure 19).
11. Proceed in a similar manner at the opposite end of the incision, impaling the vessel wall and intimal lining upon the pins situated nearest the incision end (See Figure 20a for

COUPLER sizes 3.0mm and smaller; See Figure 20b for COUPLER sizes 3.5mm and larger). Complete the pinning procedure by everting the vessel upon the remaining pins (See Figure 21a for COUPLER sizes 3.0mm and smaller; See Figure 21b for COUPLER sizes 3.5mm and larger). Ensure that both the vessel wall and the intimal layer are fully impaled upon each pin.

12. Bring the rings together by turning the Anastomotic Instrument knob clockwise, ONLY until the ejector rod has just started to move the now joined rings. Hold the Anastomotic Instrument such that the “end” pinned ring is brought to the “side” pinned ring during the approximation (See Figure 22).
13. **Prior to ejecting the joined rings, squeeze the end of the apposed jaws with a small hemostat** (See Figure 23) **to ensure ring approximation and a tight friction fit.** Turn the Anastomotic Instrument knob further clockwise to eject the joined rings.
14. Check the anastomosis under the operating microscope before opening the vascular clamps. Remove the clamps and **inspect the anastomotic site to ensure that the anastomosis has been satisfactorily completed** (patent vessel without leakage).
15. To remove the jaw assembly turn the Anastomotic Instrument knob fully counterclockwise (See Figure 13). Press the release button, located near the arrow on the Anastomotic Instrument, and remove the jaw assembly (See Figure 14).
16. Rinse the Anastomotic Instrument with water after use.

Special Instructions

WARNINGS

- The Anastomotic Instrument, Vessel Measuring Gauge, COUPLER Forceps, and Sterilization Tray are supplied non-sterile and must be sterilized prior to use.
- The Anastomotic Instrument, Vessel Measuring Gauge, COUPLER Forceps, and Sterilization Tray should be thoroughly inspected before use. Instruments that are damaged and/or in need of repair should not be used.

Limitations on reprocessing

No particular limitations

Instructions

Point of use:	Rinse all instruments with water after use.
Preparation for cleaning:	<ol style="list-style-type: none"> 1. Using a neutral (pH 7-10) detergent, wash each tool clean of all blood and debris after every use. Scrub each tool with a soft brush. Pay particular attention to areas where debris can accumulate. CAUTION: Use of a cleaner with a pH greater than 10 will remove the anodized layer of the Anastomotic Instrument and the Sterilization Tray. 2. Avoid use of any harsh material that can scratch or mar the surface of the instruments.

	3. Rinse the instruments thoroughly with running water. Apply a fine jet stream through the hole in the Anastomotic Instrument knob end and press the release button while rinsing to ensure that all surfaces of the instrument are cleaned.
Cleaning: Automated	Using an automated washing machine, clean at a temperature of 45° – 55° C, with a neutral (pH 7-10) cleaning solution for at least 10 minutes.
Cleaning: Manual	Place the instruments in an ultrasonic cleaner utilizing a neutral (pH 7-10) cleaning solution and clean ultrasonically for 15 minutes. Rinse the Anastomotic Instrument thoroughly again, applying a fine jet stream of water through the hole in the Anastomotic Instrument knob end.
Disinfection:	(Optional) Using an automated washing machine, thermally disinfect the Instrument at a temperature of 90° - 95° C, for a minimum of 5 minutes.
Drying:	Following cleaning by either manual or automated cleaning methods ensure instruments are fully dry. Do not exceed 100° C for 30 minutes.
Maintenance, Inspection and Testing:	<ul style="list-style-type: none"> • Ensure that all visible debris is removed to assure the continued quality of the instruments. • Lubricate the cleaned Anastomotic Instrument (including knob) with a water-soluble lubricant prior to sterilization. Failure to clean and lubricate the Anastomotic Instrument as directed may result in instrument failure.

EN

Packaging:	Package the instruments using the appropriate method for the sterilization cycle chosen.	
Sterilization:	GRAVITY AUTOCLAVES	
	Temperature	Recommended Exposure Time (not Total Cycle Time)
	270°F (132°C)	15 minutes (wrapped)
	PREVAC AUTOCLAVES	
	Temperature	Recommended Exposure Time (not Total Cycle Time)
270°F (132°C)	4 minutes (wrapped)	
	273°F (134°C)	3 minutes (wrapped)
	NOTE: It is recommended that each institution establish the efficacy of its sterilization procedure.	
Storage:	Recommended storage at controlled room temperature 20-25°C (68-77°F).	

Be certain that all delicate tipped instruments are covered with a device designed for instrument protection during storage.

SERVICE:

For Customer or Technical service, contact:
 Phone: + 205.941.0111 or 1.800.510.3318 (U.S. only)
 Fax: + 205.941.1522
 Website: synovismicro.com

Any component of the COUPLER system must be repaired, refurbished, or sharpened only by Synovis Micro Companies Alliance. Engraving or servicing of the components by an unauthorized facility will invalidate the warranty.

ACCESSORIES & PARTS:

ITEM	REF
COUPLER, 1.0MM, 1PK	GEM2750/I
COUPLER, 1.0MM, 6PK	GEM2750
COUPLER, 1.5MM, 1PK	GEM2751/I
COUPLER, 1.5MM, 6PK	GEM2751
COUPLER, 2.0MM, 1PK	GEM2752/I
COUPLER, 2.0MM, 6PK	GEM2752
COUPLER, 2.5MM, 1PK	GEM2753/I
COUPLER, 2.5MM, 6PK	GEM2753
COUPLER, 3.0MM, 1PK	GEM2754/I
COUPLER, 3.0MM, 6PK	GEM2754
COUPLER, 3.5MM, 1PK	GEM2755/I
COUPLER, 3.5MM, 6PK	GEM2755
COUPLER, 4.0MM, 1PK	GEM2756/I
COUPLER, 4.0MM, 6PK	GEM2756
MAC SYSTEM 2 C FORCEPS	GEM2741/CC
STERILIZATION TRAY, MAC	GEM2745
ANASTOMOTIC INSTRUMENT, TI	GEM2740
VESSEL MEASURING GAUGE	GEM2749
FORCEPS, COUPLER	GEM4183C

DISCLAIMER OF WARRANTY:

Synovis Micro Companies Alliance, Inc., a subsidiary of Baxter International Inc. warrants that reasonable care has been used in the manufacture of this device. This warranty is exclusive and in lieu of all other warranties whether expressed, implied, written or oral, including, but not limited to, any implied warranties of merchantability or fitness. Since SMCA has no control over the conditions under which the device is used, diagnosis of the patient, methods of administration or its handling after it leaves its possession, SMCA does not warrant either a good effect or against an ill effect following its use. The manufacturer shall not be liable for any incidental or consequential loss, damage or expense arising directly or indirectly from the use of this device. SMCA will replace any device which is defective at the time of shipment. No representative of SMCA may change any of the foregoing or assume any additional liability or responsibility in connection with this device.

REFERENCES:

1. DeLacure M and Wang H: Magnetic Resonance Imaging Assessment of a Microvascular Anastomotic Device for Ferromagnetism. *Journal of Reconstructive Microsurgery* 13:8, 1997.
2. Caution needed when performing MRI scans on patients with aneurysm clips. *FDA Medical Bulletin* Volume 23, Number 2, June 1993.
3. Synovis MCA Internal report, V-0718R rev A, Coupler MR Compatibility Verification Report.

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ARTICLE	REF
COUPLER, 1,0 MM, 1PQ	GEM2750/I
COUPLER, 1,0 MM, 6PQ	GEM2750
COUPLER, 1,5 MM, 1PQ	GEM2751/I
COUPLER, 1,5 MM, 6PQ	GEM2751
COUPLER, 2,0 MM, 1PQ	GEM2752/I
COUPLER, 2,0 MM, 6PQ	GEM2752
COUPLER, 2,5 MM, 1PQ	GEM2753/I
COUPLER, 2,5 MM, 6PQ	GEM2753
COUPLER, 3,0 MM, 1PQ	GEM2754/I
COUPLER, 3,0 MM, 6PQ	GEM2754
COUPLER, 3,5 MM, 1PQ	GEM2755/I
COUPLER, 3,5 MM, 6PQ	GEM2755
COUPLER, 4,0 MM, 1PQ	GEM2756/I
COUPLER, 4,0 MM, 6PQ	GEM2756
PINCES MAC SYSTEM 2 C	GEM2741CC
PLATEAU DE STÉRILISATION, MAC	GEM2745
INSTRUMENT ANASTOMOTIQUE, TI	GEM2740
JAUGE DE MESURE DES VAISSEAUX	GEM2749
PINCES, COUPLER	GEM4183C

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ARTIKEL	REF
COUPLER, 1,0 MM, 1-PACK	GEM2750/I
COUPLER, 1,0 MM, 6-PACK	GEM2750
COUPLER, 1,5 MM, 1-PACK	GEM2751/I
COUPLER, 1,5 MM, 6-PACK	GEM2751
COUPLER, 2,0 MM, 1-PACK	GEM2752/I
COUPLER, 2,0 MM, 6-PACK	GEM2752
COUPLER, 2,5 MM, 1-PACK	GEM2753/I
COUPLER, 2,5 MM, 6-PACK	GEM2753
COUPLER, 3,0 MM, 1-PACK	GEM2754/I
COUPLER, 3,0 MM, 6-PACK	GEM2754
COUPLER, 3,5 MM, 1-PACK	GEM2755/I
COUPLER, 3,5 MM, 6-PACK	GEM2755
COUPLER, 4,0 MM, 1-PACK	GEM2756/I
COUPLER, 4,0 MM, 6-PACK	GEM2756
MAC SYSTEM 2 C ZANGE	GEM2741CC
STERILISATIONSBEHÄLTER, MAC	GEM2745
ANASTOMOSEINSTRUMENT, TI	GEM2740
MESSINSTRUMENT FÜR BLUTGEFÄSSE	GEM2749
ZANGE, COUPLER	GEM4183C

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1. DeLacure M and Wang H: Magnetic Resonance Imaging Assessment of a Microvascular Anastomotic Device for Ferromagnetism. Journal of Reconstructive Microsurgery 13:8, 1997.
2. Caution needed when performing MRI scans on patients with aneurysm clips. FDA Medical Bulletin Volume 23, Number 2, June 1993.
3. Synovis MCA Internal report, V-0718R rev A, Coupler MR Compatibility Verification Report

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1. DeLacure M and Wang H: Magnetic Resonance Imaging Assessment of a Microvascular Anastomotic Device for Ferromagnetism. Journal of Reconstructive Microsurgery 13:8, 1997.
2. Caution needed when performing MRI scans on patients with aneurysm clips. FDA Medical Bulletin Volume 23, Number 2, June 1993.
3. Synovis MCA Internal report, V-0718R rev A, Coupler MR Compatibility Verification Report.

This page intentionally left blank.

Baxter Healthcare SA
Thurgauerstrasse 130
8152 Glattpark (Opfikon)
Switzerland

 Synovis Micro Companies Alliance, Inc. (a subsidiary of Baxter International Inc.),
439 Industrial Lane, Birmingham, AL 35211-4464 USA
(Tel) **205.941.0111** (Toll free) **800.510.3318** (Fax) **205.941.1522** www.synovismicro.com

The GEM design is a trademark of Baxter International Inc., or its subsidiaries.

0742112A
2020-05-29