
GEM

Global Excellence in Microsurgery

en	GEM BIOVER Microvascular Clamps	Instructions for use.....	2
fr	Clamps microvasculaires GEM BIOVER	Mode d'emploi	6
de	GEM BIOVER mikrovaskuläre Klemmen	Gebrauchsanweisung	10
es	Pinzas microvasculares GEM BIOVER	Instrucciones de uso	14
it	Clamp microvascolari GEM BIOVER	Istruzioni per l'uso	18
pt	Pinças Microvasculares GEM BIOVER	Instruções de utilização	22
nl	GEM BIOVER Microvasculaire klemmen	Gebruiksaanwijzing.....	26
pl	Zaciski mikronaczyniowe GEM BIOVER	Instrukcja użycia.....	30
sv	GEM BIOVER mikrovaskulära klämmor	Användarinstruktioner	34
no	GEM BIOVER mikrovaskulære klemmer	Instruksjoner for bruk	38
da	GEM BIOVER Mikrovasculære klemmer	Betjeningsvejledning	42
fi	GEM BIOVER -mikrovaskulaaripuristimet	Käyttöohjeet.....	46
lv	GEM BIOVER mikrovaskulārās spaiļes	Lietošanas instrukcija.....	50
zh-CN	GEM BIOVER 微血管钳	使用说明	54
zh-TW	GEM BIOVER 微血管鉗	使用說明	58

SYMBOL GLOSSARY

Symbol Glossary per US FD&C Act:

Standard	Symbol	Symbol Title	Symbol Meaning	Symbol Number
ISO-15223-1*		Do not re-use	Do not re-use	5.4.2
ISO-15223-1		Sterilized using irradiation	Sterilized using irradiation	5.2.4
ISO-15223-1		Do not use if package is damaged	Do not use if the product sterile barrier or its packaging is compromised	5.2.8
ISO-15223-1		Double sterile barrier system	Indicates two sterile barrier system	5.2.12
ISO-15223-1		Consult instructions for use	Consult instructions for use	5.4.3
			Content	
			CAUTION: Federal (USA) Law restricts this device to sale by or on the order of a physician.	
ISO-15223-1		Catalog number	Catalog number	5.1.6
ISO-15223-1		Manufacturer	Manufacturer	5.1.1

Standard	Symbol	Symbol Title	Symbol Meaning	Symbol Number
ISO-15223-1		Do not re-sterilize	Do not re-sterilize	5.2.6
ISO-15223-1		Batch code	Lot number	5.1.5
ISO-15223-1		Use by Date	Use by Date	5.1.4
ISO-15223-1		Country of Manufacture	Identify the country of manufacture of product. The date of manufacture is adjacent to this symbol.	5.1.11
ISO-15223-1		Medical Device	Indicates the item is a medical device	5.7.7
ISO-15223-1		Unique Device Identifier	Contains unique device identifier information	5.7.10

*ISO 15223-1 Medical Devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements

Additional symbols and graphics on the product labeling that are not required by the US FD&C Act:

Graphic	Graphic Meaning	Graphic	Graphic Meaning
	Diameter		Clamping Force
	Microvascular Clamp (where indicated with color, green clamp for vein, yellow clamp for artery)		Internal tracking number
	Manufacturer part number		See IFU for symbol definition

DESCRIPTION

GEM BIOVER disposable microvascular clamps can only be used **ONE TIME**.

Two types of microvascular clamps are available:

Clamps for arteries are yellow (single (1) and double (2))		
Vessel Diameter Ø	REF	REF
<1 mm	single TKS-1	double TKS-2
1 mm – 2 mm	single TKM-1	double TKM-2
2 mm – 4 mm	single TKL-1	double TKL-2

Clamps for veins are green (single (1) and double (2))		
Vessel Diameter Ø	REF	REF
0.2 – 1 mm	single TKFV-1	double TKFV-2
<1 mm	single TKSV-1	double TKSV-2
1 mm – 2 mm	single TKMV-1	double TKMV-2
2 mm	single TKLV-1	double TKLV-2

- The GEM BIOVER disposable microvascular clamps can be easily identified in the surgical site through their yellow (arteries) or green (veins) color. Yet the clamps do not reflect the light and are therefore easier on the eyes during long surgical procedures.
- GEM BIOVER disposable microvascular clamps are radiolucent, except for the axis and the spring. Both types of disposable microvascular clamps can be seen very weakly in an MRI and in fluoroscopy.
- GEM BIOVER disposable microvascular clamps deliver a precise and reliable clamping force.
- The type of vessel (artery or vein) and the vessel diameter in mm indicates to the surgeon which type of disposable microvascular clamp should be used and which clamping force is needed to be less traumatic to the vessel.

INDICATIONS FOR USE

GEM BIOVER disposable microvascular clamps are instruments which are used for all microsurgical procedures. They are used to occlude vessels during anastomosis, which is necessary as a result of vessel damage or thrombosis.

CONTRAINDICATIONS

None known.

ADVERSE REACTIONS

None known.

WARNINGS

- **Never use damaged GEM BIOVER disposable microvascular clamps and/or damaged packaging.** Use of damaged product can result in trauma of the blood vessel.
- **Use the GEM BIOVER disposable microvascular clamps only ONE TIME, SINGLE USE.** The GEM BIOVER clamps are intended for single-use only. Re-use of clamps can result in minor or major infection.
- **Do not leave any GEM BIOVER disposable microvascular clamps in the wound after finishing the surgical procedures.** Leaving clamps in the wound after surgical procedure can result in minor or major infection or trauma to tissue or vessels.
- **Do not use the GEM BIOVER disposable microvascular clamps for the central blood circulation and do not use for the central nervous system.** GEM BIOVER clamps are not validated for efficacy or endotoxin levels. The GEM BIOVER clamps are not intended for use in the central blood circulatory or central nervous system.

CAUTIONS

- **Before using any type of GEM BIOVER disposable microvascular clamps, identify the type of vessel and vessel size to make sure that the right clamp with the correct clamping force is used.** Using the wrong size or

type of clamp can result in trauma to the blood vessel.


- **The GEM BIOVER disposable microvascular clamps are gamma sterilized and cannot be resterilized. The manufacturer does not accept any liability for clamps which are resterilized by the user.** The clamps are intended for single-use only. Attempting to resterilize the clamp can result in damage to the product, trauma to the blood vessel, non-sterile product and infection.

INSTRUCTIONS FOR USE

- Choose between a single and a double GEM BIOVER microvascular clamp.
- Use only products which are in the original and undamaged packaging to ensure that the sterility is guaranteed.
- Check the expiry date on the package to make sure that the product is still sterile.
- Place blood vessel approximately 1/3 of the way into the jaws of the clamp for ideal clamping force.
- Before closing the surgical site, check that all clamps have been removed.
- After any surgical procedure, the empty sterile pouches must be counted in order to determine the precise number of GEM BIOVER disposable microvascular clamps used.
- Dispose of used device per local ordinances.

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Biover and the GEM design are trademarks of Baxter International Inc., or its subsidiaries.
Synovis Micro Companies Alliance, Inc. (a subsidiary of Baxter International Inc.)