

Biover

Biover Disposable Microvascular Clamps Consistent Pressure Every Time.

Your choice of a vessel clamp should result in secure, consistent, low-pressure occlusion. Slipping and vessel trauma from uneven or excessive clamping forces can cause complications during a procedure.

Our Biover Microvascular Clamps are unique. Their advanced hemispheric design and high-quality polycarbonate gripping surfaces ensure a secure non-slip occlusion and deliver precise, uniform, atraumatic pressure.

Upgrade your microsurgical toolkit with Biover Microvascular Clamps, from Synovis Micro Companies Alliance, Inc. – as reliably precise and consistent as world-class surgery itself.

Ordering Information

Arterial Clamps

- TKS-1** Small, single clamp - 40g/mm² clamping force
- TKS-2** Small, double clamp - 40g/mm² clamping force
- TKM-1** Medium, single clamp - 60g/mm² clamping force
- TKM-2** Medium, double clamp - 60g/mm² clamping force
- TKL-1** Large, single clamp - 120g/mm² clamping force
- TKL-2** Large, double clamp - 120g/mm² clamping force

Venous Clamps

- TKFV-1** Micro, single clamp - 15g/mm² clamping force
- TKFV-2** Micro, double clamp - 15g/mm² clamping force
- TKSV-1** Small, single clamp - 20g/mm² clamping force
- TKSV-2** Small, double clamp - 20g/mm² clamping force
- TKMV-1** Medium, single clamp - 30g/mm² clamping force
- TKMV-2** Medium, double clamp - 30g/mm² clamping force
- TKLV-1** Large, single clamp - 40g/mm² clamping force
- TKLV-2** Large, double clamp - 40g/mm² clamping force

Rx Only. See instructions for use for indications, warnings, precautions, side effects and contraindications.
Biover Microvascular Clamps are distributed in the USA by Synovis Micro Companies Alliance, Inc.

