

GEM SuperFine MicroClip

 **GEM Precision Tools**
for the enhancement of microsurgery

Another little GEM just got smaller!

The GEM titanium, hemostatic SuperFine MicroClip (GEM) is over 35% smaller than the GEM MicroClip. The secure closure of the SuperFine MicroClip is formed by its Chevron shape. The tip-to-tip closure is progressive and encompassing, and the atraumatic triangular design of the SuperFine MicroClip ensures a uniform closure. The inside of the SuperFine MicroClips are lined with diamond-shaped grooving that create interlocking teeth that enhance security and vessel closure while preventing multi-directional slippage. GEM SuperFine MicroClips are great for side branch occlusions.

Features of the GEM SuperFine MicroClip

- GEM SuperFine MicroClip cartridge ensures easy loading
- GEM SuperFine MicroClips have an optimal, atraumatic clip design
- GEM SuperFine MicroClips are pure, strong, malleable titanium
- GEM SuperFine MicroClips are MRI compatible


The GEM SuperFine MicroClip System consists of:

GEM1521 SuperFine MicroClips; 30 blue cartridges per box

GEM715 SuperFine MicroClip Applier; 15 cm length (silver ring handles)

GEM719 SuperFine MicroClip Applier; 19 cm length (silver ring handles)

Please visit our website for additional products at www.synovismicro.com



The GEM SuperFine MicroClip is manufactured by Vitalitec International and distributed by Synovis Micro Companies Alliance, Inc.

To place an order please contact Synovis MCA at 1-800-510-3318 or orders@synovismicro.com

Synovis[®]
Micro Companies Alliance, Inc.

"The Microsurgeon's Most Trusted Resource[®]"