

MERCIAN

Surgical Instruments of Excellence

Today's Microsurgeon Brings an Outstanding Background to the Table... And So Do We.

Mercian Visibility Background Material provides improved visualization of delicate nerves and vessels during the most critical part of your case. Better visualization means less fatigue and less frustration.

One side of these reversible 50mm x 25mm silicone strips – blue, green or yellow – displays a 1mm grid, for a sense of scale. The reverse side is plain.


The strips can be trimmed to accommodate almost any surgical site. The non-reflective surface eliminates distractions from theater lighting, and the non-tacky finish aids in mobilizing the most delicate nerves and vessels.

Radiopaque and latex-free, the Mercian Visibility Background Material from Synovis Micro Companies Alliance is an outstanding improvement for the microsurgeon.

Mercian Visibility Background Mats – While you'll see the advantages immediately, we still invite you to take a look at our background under the microscope!

Ordering Information

- VB2** Yellow (Box of 12)
- VB3** Green (Box of 12)
- VB4** Blue (Box of 12)



Mercian Visibility Background Material is distributed in the USA by Synovis Micro Companies Alliance, Inc.

42002 09/10